SARS-CoV-2 Antigen Rapid Detection Kit
Instructions for Use

INTENDED USE

CAS-Envision SARS-CoV-2 Antigen Rapid Detection Kit is a colloidal gold immunochromatography intended for the qualitative detection of nucleocapsid antigens from SARS-CoV-2 in human nasal vestibule swab or oral cavity and saliva swab specimen. The test provides preliminary test results but cannot be used as the sole basis for treatment or other management decision.

SUMMARY

Coronavirus is a large virus family, known to cause colds and more serious diseases such as Middle East Respiratory Syndrome (MERS) and Severe Acute Respiratory Syndrome (SARS). The new coronavirus is a new strain of coronavirus that has never been found in humans before.

The main transmission routes of SARS-CoV-2 are respiratory droplets and contact transmission. The transmission routes of Relative sensitivity aerosol and feces-oral still need to be further clarified.

PRINCIPLE

CAS-Envision SARS-CoV-2 Antigen Rapid Detection Kit (Lateral Flow Method) is based on the principle of Immunochromatography sandwich for determination of SARS-CoV-2 antigen extracted from the nasal vestibule swab or oral cavity and saliva swab specimen.

indicating that the test has been performed correctly. When the SARS-CoV-2 antigen level in the specimen is zero or below the target threshold, the T line is not colored, this indicates a negative result.

WARNINGS AND PRECAUTIONS

- 1. This kit is for in vitro diagnostic use only;
2. Please follow the correct process for specimen collection and testing to avoid false or invalid results;
3. Do not reuse any kit components;
4. Please check the packaging before use. If it is damaged, poorly sealed or has expired, please do not use it;

MATERIALS PROVIDED

Table with 4 columns: Components, Kit size, 1T/kits, 5T/kits, 25T/kits. Rows include Kit size (Foil bag, Extraction tube, Sterile swab, Biohazard bag, Instructions, QC card) and Components.

STORAGE AND STABILITY

- 1. Store at 2~30°C;
2. The test cartridge should be used within 1 hour after taking out from the foiled pouch;
3. Keep away from sunlight, moisture and heat;
4. The product batch number, production and expiry date are printed on the outer packaging box.

CAUTION BEFORE TEST

- 1. The test should be operated at room temperature (15-30 °C).
2. Please read the instructions for use carefully before performing the test.
3. Do not read results after 30 minutes.

SPECIMEN COLLECTION


Specimens collection may significantly affect the test result. Specimens should be collected carefully follow the following methods or infection control

procedures. Choose one of the following two ways to collect specimens as required.
Nasal vestibule swab specimen collection:

- 1. Take out the sample collect swab included with the kit and carefully insert it into your nostrils 2-4cm until resistance is encountered;
2. Roll the swab against the nasal mucosa and wipe it back and forth 5 times to ensure adequate sample is collected;
3. Using the same swab to collect the sample from the other nostril with the same way;

- Oral cavity and saliva swab specimen collection:
1. Do not put anything in your mouth, including food, beverages, medicine, chewing gum, and tobacco products, for at least 30 minutes before the sample is collected;
2. Cough twice before collecting the sample;

Before the test, use the browser to scan the QR code on the front of the outer box to watch the demo video and more detailed information


SARS-CoV-2 Antigen Rapid Detection Kit

TEST PROCEDURE

- 1. Take out the extraction tubes(with prefilled buffer)of the tested quantity and place them on a suitable test tube rack with the opening facing up, tear off the aluminum film;
2. Put the collected sample swab into the extraction tube, press the swab head firmly on the tube wall, and rotate the swab at least 5 times to release the antigen from the swab head into the extraction solution;
3. Keep down the swab in the extraction tube and let it stand for 1 minute;
4. Remove the swab head from the extraction solution to close to the edge of the extraction tube, and squeeze the swab head as much as possible to dry it, to let more sample remain in the extraction tube;

RESULT INTERPRETATION

- Positive Result: Colored bands appear at both test line (T) and control line (C) showing as following picture. It indicates a positive result for the SARS-CoV-2 antigen in the specimen.
Negative Result: Colored band appear at control line (C) only. It indicates that the concentration of the SARS-CoV-2 antigen is zero or below the detection limit of the test.
Invalid Result: No visible colored band appear at control line after performing the test. The directions may not have been followed correctly or the test may have deteriorated.

REPORTING OF RESULTS

Positive Test: Positive for the presence of SARS-CoV-2 antigen. Positive results indicate the presence of viral antigens, you have a very high probability of contracting SARS-CoV-2.
Negative Test: Negative results are presumptive. You have a low probability of contracting SARS-CoV-2. But negative test results do not preclude infection.

Invalid Test: The invalid result indicates that this test is invalid. Please read the instruction carefully again and repeat the test. Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control band failure.

QUALITY CONTROL

A procedural control is included in the test. A colored line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient liquid volume, adequate membrane wicking and correct procedural technique.

PERFORMANCE VERIFICATION

Sensitivity and Specificity: The nasal swabs (anterior nasal cavity) and oral saliva swabs of 118 patients with novel Coronavirus nucleic acid positive and 127 subjects with novel Coronavirus nucleic acid negative were tested simultaneously using the product and PCR method to evaluate the sensitivity, specificity and accuracy of the product;
The nasal and oral swabs test results were compared with PCR, the sensitivity, specificity and accuracy were calculated.

Table showing Relative sensitivity (97.5%), Relative specificity (99.2%), and Accuracy (98.4%) for Nasal swabs test data statistics.

Table showing Relative sensitivity (95.8%), Relative specificity (98.4%), and Accuracy (97.1%) for Oral saliva swabs test data statistics.

Minimum detection limit: Comparison test with The PCR kits, the LOD range finding was located as: 1:10000 titer of the contrived sample(the Ct value = 34), which the three repeated tests result was 100% +(F). When 1:100000 titer of the contrived sample, the Ct value of 36.227 at which the result was +(F) and 67.7%. So ct value of 34 was confirmed as the LOD in nasal swabs sample.

Cross-reactivity: Professional test, there is no cross-reaction between the pathogen samples listed in the table below and SARS-CoV-2 Antigen Rapid Detection Kit.

Table with 4 columns: S.N., Potential Cross Reactant, Concentration Tested, Cross Reactivity (Yse/ No). Lists various pathogens like Human coronavirus 299E, Human coronavirus OC43, etc.

SARS-CoV-2 Antigen Rapid Detection Kit

HOOK EFFECT

Within the titer range of clinically positive samples of SARS-CoV-2 antigens, there is no hook effect in the test results of this product.

INDEX OF SYMBOL

Table mapping symbols to meanings: IVD (In Vitro Diagnostic Use), Tests per Kit, Batch Number, Manufacturer, Expiry Date, Keep Dry, Store between 2~30°C, See Instruction for Use, Manufacturing Date, Authorized Representative, Do not reuse, Keep away from Sunlight, Catalog #, This Side Up.

LIMITATIONS OF PROCEDURE

- 1. This product can only be used to detect the N antigen of the SARS-CoV-2 in human nasal vestibule swab or oral cavity and saliva swab.
2. Positive test results do not rule out co-infections with other pathogens.
3. Sensitivity of the test after the first five days of the onset of symptoms has been demonstrated to decrease as compared to a RT-qPCR SARS-CoV-2 assay.

GENERAL INFORMATION

Shenzhen CAS-Envision Medical Technology Co., Ltd. Add: 4F, Bldg 4, Tusincere Technology Park, NO.333 Longfei Avenue, Longcheng, Longgang, Shenzhen, Guangdong, China. Website: www.cas-envision.com

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